



Instructions for Collection:

- Do not brush teeth before collecting sample.
- Do not eat or drink for 15 minutes prior to sample collection.
- Collect sample between 11:00 pm and midnight. No earlier than 11:00 pm and no later than 12:00 am, midnight, otherwise you will have to repeat the test.

This is a time sensitive collection.

It must be done between 11:00 pm and 12:00 am, midnight.

To use the Salivette:

- Remove top cap of the tube to expose the swab.
- Pour swab directly into the mouth by tipping the tube so the swab falls into the mouth.
- Do NOT touch the swab with your fingers.
- Keep the swab in your mouth for approximately one minute.
- Roll the swab in your mouth.
- Do NOT chew the swab.
- Spit the swab back into the tube.
- Do NOT touch the swab with your fingers.
- Replace the cap.
- Make sure cap is screwed on tightly.
- Record the collection time on the tube.

Additional Instructions:

If you are not able to return within eight hours after collection, please keep refrigerated until you can bring back; must be returned within three days.

You must have a lab appointment scheduled to bring in sample.
Please call two days in advance to schedule this appointment.